




Ceramics










Recommended machines and additional consumables (not included)

 CUTTING	Equipment ATM Brillant	Consumables Cut-off wheel: diamond, metal bond (bronze) Anti-corrosion coolant
 MOUNTING	Equipment Vacuum	Consumables Cold mounting: KEM 90 (porous material) KEM 35 (high density, solid material)
 GRINDING/ POLISHING	Sample size Ø 40 mm	

Pressure parameters and specimen size





Specimen diameter [mm]	25	30	40	50	60
Divergence in pressure used in the preparation methods	-(5 N...10 N)	-5 N	0	+5 N	+(5 N...10 N)

Notes:

STEP	MEDIUM		 rpm		 Single Pressure N	 min
 Planar grinding	GALAXY red	H ₂ O	250-300	▶▶ Synchronous Rotation	30	Until plane
 Grinding	GALAXY blue	H ₂ O	250-300	▶▶ Synchronous Rotation	30	2:00
 Pre-polishing	BETA	Dia-Complete Poly, 9 µm + diamond paste	120-150	◀◀ Counter Rotation	40	10:00
 Final polishing	GAMMA	Dia-Complete Poly, 3 µm + diamond paste*	120-150	▶▶ Synchronous Rotation	30	20:00

*Time and force might be divergent according to the sample size, disc size, material

BEGINNERS GUIDE

 CUTTING	<ul style="list-style-type: none"> Use suitable cut-off wheels for ceramics (e.g. diamond; metal bond) Constant cutting speed max. 0.25 mm/s
 MOUNTING	<ul style="list-style-type: none"> Use mounting material with high edge retention Cold mounting
 GRINDING	<ul style="list-style-type: none"> Start grinding with GALAXY red Continue with GALAXY blue Thoroughly wash samples and holder under running water after each grinding step
 POLISHING	<ul style="list-style-type: none"> Rinse the polishing discs with water and spin dry after use Do not stack discs with different diamond sizes Clean samples, holders and hands under running water before each polishing step Use ethanol and blow dryer to avoid water stains Check after each step under the microscope if polishing marks are of equal size and randomly oriented Rinse the OMEGA disc with water and spin dry after use Use the consumables only for ceramics and not for other materials

Notes:

SAMPLE MICROGRAPHS

OK Sample polished

20x micrograph of ceramic after GAMMA polishing

- No traces of scratches
- Clear structure/contour of the different phases

NOK Sample polished

20x micrograph of ceramic after GAMMA polishing

- Surface is not enough polished
» Repeat GAMMA step as long as necessary

Notes:

