

Copper and copper alloys

Recommended machines and additional consumables (not included)

CUTTING	Equipment ATM Brillant	Consumables Cut-off wheel: corundum, resin bond Anti-corrosion coolant
MOUNTING	Equipment ATM Opal	Consumables Hot mounting: Bakelite red/black, Thermoplast Cold mounting: KEM 20, KEM 30 Hot mounting preferred
GRINDING/ POLISHING	Sample size Ø 40 mm	

Pressure parameters and specimen size

Specimen diameter [mm]	25	30	40	50	60
Divergence in pressure used in the preparation methods	-(5 N...10 N)	-5 N	0	+5 N	+(5 N...10 N)

Notes:

STEP	MEDIUM		rpm		Single Pressure N	min
Planar grinding	SiC-paper/foil P320 (280)	H ₂ O	250-300	▶▶ Synchronous Rotation	30	Until plane
Pre-polishing	BETA	Dia-Complete Poly, 9 µm	120-150	▶▶ Synchronous Rotation	30	3:00-4:00
Polishing	SIGMA	Dia-Complete Poly, 3 µm	120-150	▶▶ Synchronous Rotation	30	3:00-4:00
Final polishing	OMEGA	Eposil F 0.1 µm*	120-150	◀◀ Counter Rotation	15	1:00-2:00* (H ₂ O during final 0:30)
Optional etching (chem.)	Cu Etchant A (chloride version)**					Approx. 0:02

* 50 ml Eposil F + 1 ml H₂O₂ + 1 ml NH₃, otherwise polishing time x2

** ATM Item No. 95000508

BEGINNERS GUIDE

CUTTING	<ul style="list-style-type: none"> Use suitable cut-off wheels (e.g. ATM NF-A wheels) Constant cutting speed max. 0.25 mm/s
MOUNTING	<ul style="list-style-type: none"> Use mounting material with high edge retention
GRINDING	<ul style="list-style-type: none"> Grind with SiC-paper/foil P320 (280) Thoroughly wash samples and holder under running water after each grinding step Use 1 sheet of SiC-paper/foil for maximum 4 samples
POLISHING	<ul style="list-style-type: none"> For new materials start with longest recommended step times and optimize to shorter times Soft high purity copper requires only 1/2 contact pressure (15 N) and prolonged OMEGA step (5 min) Rinse the polishing discs with water and spin dry after use Do not stack discs with different diamond sizes Clean samples, holders and hands under running water before each polishing step Use ethanol and blow dryer to avoid water stains Check after each step under the microscope if polishing marks are of equal size and randomly oriented Rinse the OMEGA disc with water and spin dry after use Use the consumables only for copper and not for other materials Rinse the cap of the Eposil F bottle after use, put cap back on


Notes:

SAMPLE MICROGRAPHS

OK Sample polished

20x micrograph of cast bronze after OMEGA polishing


- CuSn intermetallic defined
- Pores from casting visible, minimal residual scratches
- No pittings or relief from over-polishing



NOK Sample polished

10x micrograph of cast bronze after OMEGA polishing


- Scratches from 3 µm step after OMEGA
 - Clean all polishing discs with clean brush under running water
 - Repeat OMEGA step
 - Clean sample and sample holder



OK Etched Sample A

10x micrograph of cast bronze, etched with Cu Etchant A


- Good brightness contrast between regions of different dendrite orientation



OK Etched Sample B

20x micrograph of Cu, etched with Cu Etchant A after 5 min OMEGA polishing

- Grain boundaries and twinning discriminable
- Almost no residual scratches



Notes: