

Nitrided steel

Recommended machines and additional consumables (not included)

CUTTING	Equipment ATM Brillant	Consumables Cut-off wheel: corundum, resin bond Anti-corrosion coolant
MOUNTING	Equipment ATM Opal	Consumables Hot mounting: EPO black, EPO-Max* Cold mounting: KEM 15 plus
GRINDING/ POLISHING	Sample size Ø 40 mm	Hot or cold mounting

Pressure parameters and specimen size

Specimen diameter [mm]	25	30	40	50	60
Divergence in pressure used in the preparation methods	-(5 N...10 N)	-5 N	0	+5 N	+(5 N...10 N)

Notes:

STEP	MEDIUM	H ₂ O	rpm		Single Pressure N	min
Planar grinding	SiC-Paper/foil P320 (280)	H ₂ O	250-300	▶▶ Synchronous Rotation	30	Until plane
Pre-polishing	BETA	Dia-Complete Poly, 9 µm	120-150	◀◀ Counter Rotation	30	5:00
Polishing	GAMMA	Dia-Complete Poly, 3 µm	120-150	▶▶ Synchronous Rotation	30	6:00
Polishing	ZETA	Dia-Complete Poly, 1 µm	120-150	▶▶ Synchronous Rotation	30	3:00
Final polishing	OMEGA	Eposal 0.06 µm	120-150	◀◀ Counter Rotation	15	2:00 (H ₂ O during final 0:30)
Optional: Etching (chem.)	Kalling II**					Approx. 0:02-0:10

* Before hot mounting the specimen should be wrapped in aluminum foil to stabilize the nitrided layer

** ATM Item No. 92002347

BEGINNERS GUIDE

CUTTING	<ul style="list-style-type: none"> Use suitable cut-off wheels for ferrous material (e.g. ATM FS-E wheels) Constant cutting speed max. 0.25 mm/s
MOUNTING	<ul style="list-style-type: none"> Use mounting material with high edge retention Hot and cold mounting both possible
GRINDING	<ul style="list-style-type: none"> Grind with SiC-Paper/foil P320 (280) Thoroughly wash samples and holder under running water after each grinding step
POLISHING	<ul style="list-style-type: none"> Rinse the polishing discs with water and spin dry after use Do not stack discs with different diamond sizes Clean samples, holders and hands under running water before each polishing step Use ethanol and blow dryer to avoid water stains Check after each step under the microscope if polishing marks are of equal size and randomly oriented Rinse the OMEGA disc with water and spin dry after use Use the consumables only for nitrided steel and not for other materials Rinse the cap of the Eposal bottle after use, put cap back on Use cosmetic tissues to clean possible traces of Eposal after the last polishing step

Notes:

SAMPLE MICROGRAPHS

OK Sample polished

20x micrograph of nitrided steel after OMEGA polishing

- No traces of scratches
- Clear structure/contour of the different phases



NOK Sample polished

20x micrograph of nitrided steel after OMEGA polishing

- Spalling in the nitride layer (created by too much pressure)
 - » Check the chosen parameters of the final polishing step
 - » Repeat all (!) steps carefully



OK Sample etched

20x micrograph of nitrided steel etched with Nital 3%

- Nitride layer clearly visible



Notes: